




# Biospecimen Collection, Processing, and Shipment Manual

## Appendix B: IUGB Stat Lab Processing Sheet

<b>OxICU K-76</b>  Protocol: MMGE-NIA-NCRAD-OxICUK76 ARM: CASE: FORMAT: KIT NUMBER		<b>IUGB STAT LAB PROCESSING SHEET</b>				COLLECTION DATE:	
		VERSION: 5 INITIALS: DGP	LAST UPDATED: 04/22/2025	VISIT:	KIT#:	COLLECTION TIME:	
		EFFECTIVE DATE:				TIME REC'D:	
<b>FINISH ALL PROCESSING WITHIN TWO (2) HOURS OF DRAW TIME UNLESS OTHERWISE NOTED.</b>							
TUBE TYPE/ QUANTITY	BARCODES	CENTRIFUGE	HEMOLOGY /TURBIDITY SCALE 0-3	INSTRUCTIONS/ # SAMPLES	ACTUAL # ALIQUOTS (μL)	TIME FROZEN	
1 x 10 mL EDTA (Purple)	_____ BARCODE	TEMP: 4°C SPEED: 2000 x G TIME: 10 MINUTES  CENTRIFUGE START TIME:		1) VOLUME: 5 x 100uL 2 x 500uL REMAINDER AT 1500 uL KEEP RESIDUAL 2) TUBES: SARSTEDT 3) CAPS: PURPLE 4) STORE: 1 FULL ALIQUOT INTO BU, REST INTO MAIN	_____ x 100 uL ALIQUOTS	DRY ICE/-80°C	
	_____ CHILD ALIQUOTS				_____ x 500 uL ALIQUOTS		
	_____ BUFFY COAT				_____ x 1500 uL ALIQUOTS  1 x _____ uL RESIDUAL		
1 x 10 mL PAXgene	_____ Barcode			1) RECORD VOLUME 2) STORE: PAX MAIN	_____ uL	DRY ICE/-80°C	
1 x 10 mL NaHep (Green)	_____ Barcode			1) RECORD VOLUME 2) ASSIGN STORAGE: RI-IUGB AMBIENT TRANSFER FROM RILEY TO TK. ADD TO AMBIENT COC. TRANSFER TO TK VIA PATH LAB TUBE.	_____ uL		
NOTES: Day 1 – only RNA collected Day 4 – only 1 10mL EDTA collected Other Visits will have 10mL EDTA and 10 mL NaHep collected  Processed in BOTH Riley and Noyes Kit is Pre-built							
PROCESSED BY / DATE		ENTERED DATA / DATE:		LABELED BY / DATE:		STORED / DATE:	
QC'D BY / DATE:							